

UPCERA



A7+

Porcelain Furnace

Reliable Materials Expert



A7+

Porcelain Furnace

- Double-layer metal protection to the chamber
- Integrated vacuum technology
- Automatic platform lift mode
- Built-in error warning
- 2-Stage program to meet special sintering demands
- Temperature compensation
- Vacuum testing, cleaning, drying, and other maintenance procedures
- Unique, homogeneous heat distribution

Device Parameters

Specification	Parameters
Dimensions	330*420*560mm
Weight	32.5kg
Voltage	220V±10%/50Hz
Protection Level	IPX1
Protective Tube1	3.0A
Protective Tube2	8.0A
Heating Rate	5~100 C
Max Working Temperature	1200 C
Constant Temperature Time	00:30-30:00min
Chamber Available Sizes	φ85×55 (mm)
Ultimate Vacuum Degree	<35mmHg